

Xrf & Kbr Plastic Ring Lab Powder Pellet Pressing Mold For Ftir

Item Number: PMXP



Introduction

Get precise XRF samples with our plastic ring lab powder pellet pressing mold. Fast tableting speed and customizable sizes for perfect molding every time.

[Learn More](#)

Instrument model	PMXP
Press the shape of the sample	
Heating temperature	Room temperature-300C
Mold material	Alloy tool steel
Sample size	Φ25mm (d)
Sample thickness	15,25,50,100,250,500um (6 quantitative rings)
External dimension	200*60mm (D*H)
Power supply	220V/300W
Size diagram	