

## Xrf & Kbr Plastic Ring Lab Powder Pellet Pressing Mold

Item Number: PMXP



### Introduction

Get precise XRF samples with our plastic ring lab powder pellet pressing mold. Fast tableting speed and customizable sizes for perfect molding every time.

[Learn More](#)

|                     |             |
|---------------------|-------------|
| Model               | PMXP-A      |
| Material            | Cr12MoV     |
| Indenter hardness   | HRC60-HRC62 |
| Sample size         | <φ50 mm     |
| Sample thickness    | 0.02-0.1 mm |
| External dimensions | φ200×50 mm  |
| Weight(Kg)          | 3           |

**Other sizes can be customized**