

Xrf & Kbr Steel Ring Lab Powder Pellet Pressing Mold For Ftir

Item Number: PMXS



Introduction

Produce perfect XRF samples with our steel ring lab powder pellet pressing mold. Fast tableting speed and customizable sizes for accurate molding every time.

[Learn More](#)

| Instrument model | PMXS |
|--|---------------------------|
| Sample shape | |
| Die material | Alloy tool steel :Cr12MoV |
| Indenter hardness | HRC60-HRC62 |
| Sample size | Φ32[]Φ40mm (M) |
| Cavity depth | 45m (N) |
| Dimensions | Φ73*133mm(L*H) |
| Weight | 3.2Kg |
| Diagram of hydraulic powder press size | |